22_Denoising microscope images

https://www.youtube.com/watch?v=2PrzKWkqOtU&t=4s (https://www.youtube.com/watch?v=2PrzKWkqOtU&t=4s)

'Non local means' denoising filter is on scikit-image

But let's first of all, read an image

In [1]:

```
from skimage import io, img as float
```

I will read and then convert the image into float because later I will have to do some Math

In [2]:

```
img = img_as_float(io.imread("python_for_microscopists-master/images/denoising/noisy_img.jp
```


I will import pyplot to save and display images

In [3]:

```
from matplotlib import pyplot as plt
```

In [4]:

plt.imshow(img);



Now I will apply gaussian denoising from scipy

In [5]:

from scipy import ndimage as nd

In [6]:


```
gaussian_img = nd.gaussian_filter(img, sigma=3)
```

In [7]:

```
plt.imsave("python_for_microscopists-master/images/denoising/gaussianMINE.jpg", gaussian_im
```

In [8]:

plt.imshow(gaussian_img);



As you've just seen, the gaussian filter removes noise but the image becomes blurred.


So I will try now a median filter (edge preserving filter)

In [9]:

```
median_img = nd.median_filter(img, size=3)
plt.imsave("python_for_microscopists-master/images/denoising/medianMINE.jpg", median_img)
```

In [10]:

plt.imshow(median_img);



Now we want to apply the denoise filter from scikit-image.

You find all the documentation here: https://scikit-

<u>image.org/docs/dev/auto_examples/filters/plot_nonlocal_means.html (https://scikit-image.org/docs/dev/auto_examples/filters/plot_nonlocal_means.html)</u>

This is the command we are interested in:

slow algorithm

```
denoise = denoise_nl_means(noisy, h=1.15 * sigma_est, fast_mode=False, **patch_kw)
```

The denoise_nl_means filter is in skimage.restoration library

In [11]:

```
from skimage.restoration import denoise_nl_means, estimate_sigma
```

estimate sigma belongs to numpy so we have to import numpy

In [12]:

```
import numpy as np
```

In [13]:

```
# estimate the noise standard deviation from the noisy image
sigma_est = np.mean(estimate_sigma(img, multichannel=True))
```

Now I will apply the denoise_nl_means filter here:


In [14]:

In [15]:

```
plt.imsave("python_for_microscopists-master/images/denoising/nlm.jpg", nlm)
```


In [16]:

```
plt.imshow(nlm);
```



Now let's look at the original image and the denoised image side-by-side

In [18]:



In []:			